



Nitrocefin PRODUCT DATA SHEET

issue date 12/12/2018

| | |
|----------------------------|---|
| Product Name: | Nitrocefin |
| Product Number: | N005 |
| CAS Number: | 41906-86-9 |
| Molecular Formula: | $C_{21}H_{16}N_4O_8S_2$ |
| Molecular Weight: | 516.51 |
| Form: | Powder |
| Appearance: | Orange or yellow powder |
| Source: | Synthetic |
| Storage Conditions: | Protect from light. Store with inert gas. -20°C |

Description:

Nitrocefin is a cephalosporin with chromogenic properties and is routinely used to detect beta-lactamase enzymes produced by beta-lactam resistant bacteria. Nitrocefin is soluble in DMSO and is commonly used at a 1.0 mg/mL concentration.

This product is considered a dangerous good. Quantities above 1 g may be subject to additional shipping fees. Please contact us for specific questions.

Zhangming et al. used nitrocefin from TOKU-E as a substrate to study TEM-1 β -lactamase activity from *E. coli*. ["Label-Free Measurements of Reaction Kinetics Using a Droplet-Based Optofluidic Device."](#)

Liu et al. used nitrocefin from TOKU-E to study and develop a homogeneous biosensor. ["Parts-per-Million of Polyethylene Glycol as a Non-Interfering Blocking Agent for Homogeneous Biosensor Development."](#)

Dahlin et al. used nitrocefin from TOKU-E in a surrogate β -lactamase-nitrocefin assay. ["A Cell-Free Fluorometric High-Throughput Screen for Inhibitors of Rtt109-Catalyzed Histone Acetylation."](#)

Ohlhoff et al. used nitrocefin from TOKU-E as a substrate to study the activity of EstG34 β -lactamases. ["An unusual feruloyl esterase belonging to family VIII esterases and displaying a broad substrate range"](#)

Huang et al. used nitrocefin from TOKU-E as a substrate to study the activity of VIM-2 Metallo- β -lactamases (MBLs). ["Inhibiting the VIM-2 Metallo- \$\beta\$ -Lactamase by Graphene Oxide and Carbon Nanotubes."](#)

Eze E et al. used nitrocefin from TOKU-E to confirm beta-lactamase production in *E. coli* and *Klebsiella* species from Nigeria. Read more here: ["Occurrence of Beta-Lactamases and the Antibiogram Pattern of Clinical Isolates of Escherichia coli and Klebsiella Species in Nsukka Metropolis."](#)

Choi et al. measured enzymatic activity of engineered protein switches by exploiting nitrocefin hydrolysis. Read more here: ["Electrochemical Activation of Engineered Protein Switches."](#)

Pierre, et al. used nitrocefin from TOKU-E to measure the enzymatic activity of various beta-lactamases. Read more here: ["Molecular Determinants for Protein Stabilization by Insertional Fusion to a Thermophilic Host Protein."](#)

Tullman and Nicholes, et al. used nitrocefin from TOKU-E to study and characterize enzymatic protein switches. Read more here: ["Enzymatic protein switches built from paralogous input domains."](#)


Mechanism of Action:

Essentially all beta-lactamase enzymes hydrolyze the amide bond between the carbonyl carbon and the nitrogen in the beta-lactam ring of nitrocefin. Macroscopic detection of this process is made possible because a ultraviolet absorption shift from intact versus hydrolyzed nitrocefin occurs within the visible light spectrum (~380 nm to ~500 nm, or yellow to red).

Microbiology Applications

Nitrocefin is used to detect beta-lactamase activity from suspected beta-lactam resistant bacteria (see protocol below).

Technical Data:



Example of nitrocefin color change before and after exposure to beta-lactamase.

(A) Concentrated nitrocefin (10.0 mg/mL) in DMSO before dilution with PBS buffer. **(B)** Nitrocefin diluted with PBS buffer to working concentration (1.0 mg/mL). The yellow color is indicative of intact, undegraded nitrocefin. **(C)** 25 units of beta-lactamase dropped on top of nitrocefin (1.0 mg/mL in PBS). The red color is the result of beta-lactamase mediated cleavage of the nitrocefin. **(D)** Vortexed mixture of contents shown in picture (C).

References:

O'Callaghan, Cynthia H. et. al. "Novel Method for Detection of B-Lactamases by Using a Chromogenic Cephalosporin Substrate." *Antimicrobial Agents and Chemotherapy* 1.4 (1972): 283-88. *Ncbi.gov* Web. 10 Oct. 2012.

Parr, T. R., Jr. "Simple Screening Method for Beta-lactamase-positive and -negative Ampicillin-resistant Haemophilus Influenzae Isolates." *Journal of Clinical Microbiology* 20.1 (1984): 131-32. *Ncbi.gov*. Web. 10 Oct. 2012.

If you need any help, contact us: info@toku-e.com. Find more information on: www.toku-e.com/