

# 5-Fluorocytosine PRODUCT DATA SHEET

issue date 01/06/2020

**Product Name:** 5-Fluorocytosine

Product Number: F003

CAS Number: 2022-85-7 Molecular Formula:  $C_4H_4FN_3O$ 

Molecular Weight: 129.09
Form: Powder

**Appearance:** White or almost white crystalline powder

**Solubility:** Sparingly soluble in water (10.5 mg/mL). Slightly soluble in alcohol.

Source: Synthetic

**Storage Conditions:** 2-8°C; Protect from light.

**Description:** 5-Fluorocytosine (Flucytosine) is a fluorinated cytosine (pyrimidine) analog

with antifungal activity. It was first developed in 1957 and is a competitive inhibitor of uracil metabolism. 5-Fluorocytosine is taken up by the fungus and is

converted to the cytostatic 5-fluorouracil by cytosine deaminase. The

compound is active against *Candida* and *Cryptococcus*. 5-Fluorocytosine can also be used in gene therapy, using a virus to deliver the gene coding for cytosine deaminase gene to brain cancer tumors, and the resulting fluorouracil selectively kills malignant cells. 5-Fluorocytosine is sparingly soluble in water.

**Mechanism of Action:** 5-Fluorocytosine (5-FC) is taken up by the fungus and is deaminated by

cytosine deaminase (CD) and converted to the cytostatic 5-fluorouracil (5-FU). During transcription, the pyrimidine analog is incorporated into the RNA which inhibits subsequent rounds of transcription and building of certain essential

proteins.

An additional mechanism of action involves the conversion of 5-Fluorocytosine to 5-fluorodeoxyuridinemonophosphate, which inhibits fungal DNA synthesis.

**Spectrum:** 5-Fluorocytosine is active *in vitro* and *in vivo* against some strains of *Candida* 

and Cyptococcus species.

Microbiology Applications 5-Fluorocytosine was found to be very active in vitro in a study with 8803 Candida isolates representing 18 yeast species (8803 isolates) using broth microdilution testing, with 95% of the isolates being susceptibility. The compound had a prolonged post-antifungal effect and concentrationindependent activity (Pfaller et al, 2002).

> In a study of over 1000 isolates of *Candida exposed to 5-Fluorocytosine*, 83.4% of them were susceptible. Resistance was uncommon (10.4% of isolates) and an MIC > 32 µg/ml was an indication of resistance in vitro (Quindós et al. 2004).

Representative MIC data for 5-Fluorocytosine is available from our Antimicrobial Index.

## **Plant Biology Applications**

5-Fluorocytosine can be used as a negative selection agent in transgenic plants using a bacterial cytosine deaminase codA gene in tobacco plant systems. In transformed callus, expression of this gene results in cell death on 5-Fluorocytosine (Stougaard et al, 1993).

### **Cancer Applications**

5-Fluorocytosine can be used in gene therapy applications in mouse models of brain cancer. Researchers used replicating virus Toca 511 that infects dividing cells, preferentially malignant cells. The virus encodes an optimized yeast-derived cytosine deaminase (CD) that converts 5-Fluorocytosine (5-FC) into 5-fluorouracil (5-FU) which selectively kills dividing cells. This allows infected cells to remain as a reservoir of virus capable of infecting new malignant cells during subsequent viral proliferation. Authors used 2 immunocompetent mouse brain tumor models (CT26 in BALB/c mice and Tu-2449 in B6C3F1 mice) to show that the conversion proceeded in vivo. Repeated cycles of 5-FC administration was able to shrink and prevent further growth of the IC tumor (Ostertag et al, 2012).

#### References:

Marco F et al (2003) Trends in frequency and *in vitro* susceptibilities to antifungal agents, including voriconazole and anidulafungin, of *Candida* bloodstream isolates. Results from a six-year study (1996-2001) Diagn. Microbiol. Infect. Dis. 46(4):259-264 PMID 12944017

Ostertag D et al (2012) Tumor eradication and prolonged survival from intratumoral conversion of 5-fluorocytosine to 5-fluorouracil using a nonytic retroviral replicating vector. Neuro. Oncol. 14(2):145-159 PMID 22070930

Pfaller, MA et al (2002) *In vitro* activities of 5-fluorocytosine against 8,803 clinical isolates of *Candida* spp.: Global assessment of primary resistance using National Committee for Clinical Laboratory Standards susceptibility testing methods. *Antimicrob. Agents and Chemother.* 46(11):3518-3521 PMID 12384359

Quindós G (2004) *In-vitro* activity of 5-fluorocytosine against 1,021 Spanish clinical isolates of *Candida* and other medically important yeasts. Rev Iberoam. Micol. 21(2):63-69 PMID 15538829

Rice LB and Mahmoud AG (1999) Antifungal agents: Mode of action, mechanisms of resistance, and correlation of these mechanisms with bacterial resistance. Clin. Microbiol. Rev. 12(4):501-17 PMID 10515900

Stougaard J (1993) Substrate-dependent negative selection in plants using a bacterial cytosine deaminase gene. Plant J. 3(5):755-761

Te Dorsthorst, DTA et al (2002) *In vitro* interaction of Flucytosine combined with amphotericin B or fluconazole against thirty-five yeast isolates determined by both the fractional inhibitory concentration index and the response surface approach. Antimicrob. Agents Chemother. 46(9):2982-2989 PMID 12183257

Waldorf AR and Polak K (1983) Mechanisms of action of 5-Fluorocytosine. Antimicrob. Agents Chemother. 23 (1):79-85

If you need any help, contact us: info@toku-e.com. Find more information on: www.toku-e.com/