

## Sulfadiazine PRODUCT DATA SHEET

issue date 01/06/2020

Product Name: Sulfadiazine

Product Number: S012

**CAS Number:** 68-35-9

Molecular Formula:  $C_{10}H_{10}N_4O_2S$ 

Molecular Weight: 250.28 g/mol

Form: Powder

**Appearance:** White or almost white crystalline powder.

**Source:** Synthetic

Storage Conditions: 2-8°C

**Description:** Sulfadiazine is a sulfonamide antibiotic that was developed in 1940 by

Richard Roblin at the Stamford Research Laboratories of the American Cyanamid Company to be a more potent and less toxic alternative to

Sulfapyridine. Sulfadiazine is a bacteriostatic folic acid synthesis inhibitor that shows broad spectrum activity against gram-positive bacteria, gram-negative

bacteria and Chlamydia.

Sulfadiazine is a competitive inhibitor of the bacterial enzyme dihydropteroate

synthetase. This enzyme is needed for the proper processing of paraaminobenzoic acid (PABA) which is essential for folic acid synthesis. The inhibited reaction is necessary in these organisms for the synthesis of folic

acid.

Sulfadiazine is sparingly soluble in aqueous solution (0.077 mg/mL).

Synonyms: 4-Amino-N-(2-pyrimidinyl)benzenesulfonamide,  $N^1$ -(Pyrimidin-2-

yl)sulfanilamide

**Mechanism of Action:** Sulfadiazine is a competitive inhibitor of the bacterial enzyme dihydropteroate

synthetase. This enzyme is needed for the proper processing of paraaminobenzoic acid (PABA) which is essential for folic acid synthesis. The inhibited reaction is necessary in these organisms for the synthesis of folic

acid.

**Spectrum:** Sulfadiazine is a broad-spectrum antibiotic targeting a wide range of gram-

positive and gram-negative bacteria especially those which cause urinary tract

infections.

Microbiology Applications Sulfadiazine is commonly used in clinical in vitromicrobiological antimicrobial susceptibility tests (panels, discs, and MIC strips) against gram positive and gram negative microbial isolates. Medical microbiologists use AST results to recommend antibiotic treatment options for infected patients. Representative MIC values include:

- Neisseria meningitidis 0.5 μg/mL 200 μg/mL
- For a complete list of sulfadiazine MIC values, click here.

Sulfadiazine is a component of SPS agar (Sulfite Polymyxin Sulfadizine Agar), that can be used for the detection, isolation, and emuneration of Clostridium perfringens and Clostridium botulinum in food; as well as the detection of V. cholerae and V. vulnificus in the environment.

Sulfadiazine is commonly used with pyrimethamine to treat toxoplasmosis in patients with acquired immunodeficiency syndrome. It is also used to treat newborns with congenital infections. Sulfadiazine has been used to control acute infections when studying murine models of reactivated toxoplasmosis.

## **Plant Biology Applications**

Sulfadiazine is a herbicide that is able to inhibit folic acid synthesis in plants. Sulfadiazine and the resistance gene sul, create an effective selection system for the creation of transgenic tobacco and algae. Colltransformation experiments in tobacco revealed that sul is even superior to nptll, the currently most efficient selectable marker gene, and thus provides an attractive marker for the highlithroughput genetic transformation of plants and algae

## References:

"Sulfadiazine: Mechanism of Action." Cancer News Network. N.p., 2009. Web. 21 Aug. 2012.

Tabatabaei, I., Dal Bosco, C., Bednarska, M., Ruf S., Meurer, J., Bock, R. (2018) A highly efficient sulfadiazine selection system for the generation of transgenic plants and algae. Plant Biotechnol. J., https://doi.org/10.1111/pbi.13004

Minden, V., Deloy, A., Volkert, A. M., Leonhardt, S. D., & Pufal, G. (2017). Antibiotics impact plant traits, even at small concentrations. AoB PLANTS, 9(2), plx010. doi:10.1093/aobpla/plx010

If you need any help, contact us: info@toku-e.com. Find more information on: www.toku-e.com/