

Temocillin Disodium PRODUCT DATA SHEET

issue date 01/06/2020

Product Name: Temocillin Disodium

Product Number: T130

CAS Number: 61545-06-0

Molecular Formula: $C_{16}H_{18}N_2Na_2O_7S_2$

Molecular Weight: 458.42

Form: Powder

Appearance: white or almost white powder

Source: semi-synthetic

Description: Temocillin Disodium is the salt of Temocillin a semisynthetic,

β-lactamase resistant carboxypenicillin antimicrobial used as an alternative to

carbapenems that is primarily used for Gram-negative bacteria.

Mechanism of Action: Like other β-lactams, Temocillin Disodium interferes with penicillin binding

protein (PBP) activity involved in the final phase of peptidoglycan synthesis. PBB catalyzes the pentaglycine crosslink between alanine and lysine. Without a pentaglycine crosslink, the integrity of the cell wall is compromised which leads to cell lysis and death. Temocillin is resistant to most β-lactamases, but

resistance is commonly due to cells containing plasmid-encoded β-

lactamases.

Spectrum: Temocillin is active against several multi-resistant Gram-negative species

including Salmonella typhimurium, E. coli, Haemophilus influenzae and others. In contrast, Gram-positive bacteria, bacteria with altered penicillin-binding proteins, and Pseudomonas aeruginosa are not susceptible to

temocillin.

Microbiology Applications Temocillin is commonly used in clinical in vitro microbiological antimicrobial

susceptibility tests (panels, discs, and MIC strips) against Gram-negative microbial isolates. Medical microbiologists use this information to recommend

antibiotic treatment options. Representative MIC ranges include:

E. coli 1.00 μg/mL - 64.00 μg/mL

- H. influenzae β-lactamase negative: 0.25 0.50 µg/ml
- H. influenzae β-lactamase positive: 0.50 1.00 µg/ml
- S. aureus >128.00 µg/ml

For a list of MIC values, click here.

References:

Pitout JD, Sanders CC, Sanders WE (1997) Antimicrobial resistance with focus on beta-lactam resistance in Gram-negative bacilli. Am. J. Med. 103:51 PMID 9236486

Jules, K and Neu HC (1982) Antibacterial Activity and Beta-lactamase Stability of Temocillin. Antimicrob. Agents and Chemother. 22(3): 453-460

Rodriguez-Villalobos H et al (2006) *In vitro* activity of temocillin against extended spectrum beta-lactamase-producing *Escherichia coli.* J Antimicrob. Chemother. 57(4):771-774 PMID 16501056

Van Landuyt HW, Pyckavet MP, Lambert A and Boelaert J (1982) *In vitro* activity of temocillin (BRL 17421), a novel B-lactam antibiotic. Antimicrob. Agents. Chemother. 22(4):535-540 PMID 7181470

If you need any help, contact us: info@toku-e.com. Find more information on: www.toku-e.com/